



CstF-64T Monoclonal Antibody

Catalog No	YP-mAb-01632
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CSTF2T
Immunogen	The antiserum was produced against synthesized peptide derived from human CSTF2T. AA range:91-140
Specificity	CstF-64T Monoclonal Antibody detects endogenous levels of CstF-64T protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TauCstF-64;TauCstF 64;KIAA0689;CstF 64T;Cleavage stimulation factor 64 kDa subunit tau variant;CstF-64T;CSTF2T;Cleavage stimulation factor subunit 2 tau variant;CF-1 64 kDa subunit tau variant;CSTF 64 kDa subunit tau variant
Observed Band	64kD
Calculated Molecular Weight	64kD
Cell Pathway	Nucleus .
Tissue Specificity	Brain,Coronary artery,Epithelium,Testis,
Function	function:May play a significant role in AAUAAA-independent mRNA polyadenylation in germ cells. Directly involved in the binding to pre-mRNAs., PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 RRM (RNA recognition motif) domain.,
Background	function:May play a significant role in AAUAAA-independent mRNA polyadenylation in germ cells. Directly involved in the binding to pre-mRNAs., PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 RRM (RNA recognition motif) domain.,




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using CstF-64T Monoclonal Antibody