



DPF2 Monoclonal Antibody

Catalog No	YP-mAb-01668
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DPF2
Immunogen	The antiserum was produced against synthesized peptide derived from human REQU. AA range:151-200
Specificity	DPF2 Monoclonal Antibody detects endogenous levels of DPF2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BAF45D;BRG1 associated factor 45D;BRG1-associated factor 45D;Protein requiem;REQ;DPF2;Zinc finger protein ubi-d4;UBID4;Apoptosis response zinc finger protein;D4;zinc and double PHD fingers family 2
Observed Band	44kD
Calculated Molecular Weight	44kD
Cell Pathway	Nucleus . Cytoplasm .
Tissue Specificity	Ubiquitous.
Function	function:May be a transcription factor required for the apoptosis response following survival factor withdrawal from myeloid cells. Might also have a role in the development and maturation of lymphoid cells.,similarity:Belongs to the requiem/DPF family.,similarity:Contains 1 C2H2-type zinc finger.,similarity:Contains 2 PHD-type zinc fingers.,subcellular location:30% nuclear. 70% cytoplasmic.,tissue specificity:Ubiquitous.,
Background	The protein encoded by this gene is a member of the d4 domain family, characterized by a zinc finger-like structural motif. This protein functions as a transcription factor which is necessary for the apoptotic response following



deprivation of survival factors. It likely serves a regulatory role in rapid hematopoietic cell growth and turnover. This gene is considered a candidate gene for multiple endocrine neoplasia type I, an inherited cancer syndrome involving multiple parathyroid, enteropancreatic, and pituitary tumors. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using DPF2 Monoclonal Antibody