




DQX1 Monoclonal Antibody

Catalog No	YP-mAb-01669
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	DQX1
Immunogen	The antiserum was produced against synthesized peptide derived from human DQX1. AA range:551-600
Specificity	DQX1 Monoclonal Antibody detects endogenous levels of DQX1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ATP-dependent RNA helicase DQX1;DEAQ box polypeptide 1;EC:3.6.4.12; DQX1
Observed Band	79kD
Calculated Molecular Weight	79kD
Cell Pathway	Nucleus .
Tissue Specificity	Placenta,
Function	similarity:Contains 1 helicase ATP-binding domain.,similarity:Contains 1 helicase C-terminal domain.,
Background	similarity:Contains 1 helicase ATP-binding domain.,similarity:Contains 1 helicase C-terminal domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using DQX1
Monoclonal Antibody