



E4BP4 Monoclonal Antibody

Catalog No	YP-mAb-01676
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NFIL3
Immunogen	The antiserum was produced against synthesized peptide derived from human NFIL3. AA range:61-110
Specificity	E4BP4 Monoclonal Antibody detects endogenous levels of E4BP4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NFIL3A;Interleukin-3-binding protein 1;Interleukin-3 promoter transcriptional activator;IL3BP1;E4BP4;NFIL3;Nuclear factor interleukin-3-regulated protein;E4 promoter-binding protein 4;Transcriptional activator NF-IL3A
Observed Band	51kD
Calculated Molecular Weight	51kD
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in bladder stomach, thyroid, spinal cord, lymph node, trachea, adrenal gland, bone marrow and muscle.
Function	function:Acts as a transcriptional regulator that recognizes and binds to the sequence 5'-[GA]TTA[CT]GTAA[CT]-3', a sequence present in many cellular and viral promoters. Represses transcription from promoters with activating transcription factor (ATF) sites. Represses promoter activity in osteoblasts (By similarity). Represses transcriptional activity of PER1 (By similarity). Represses transcriptional activity of PER2 via the B-site on the promoter (By similarity). Activates transcription from the interleukin-3 promoter in T-cells. Competes for the same consensus-binding site with PAR DNA-binding factors (DBP, HLF and TEF) (By similarity). Component of the circadian clock that acts as a negative



regulator for the circadian expression of PER2 oscillation in the cell-autonomous core clock (By similarity). Protects pro-B cells from programmed cell death., induction:Up-regulated by PHA or TP

Background

The protein encoded by this gene is a transcriptional regulator that binds as a homodimer to activating transcription factor (ATF) sites in many cellular and viral promoters. The encoded protein represses PER1 and PER2 expression and therefore plays a role in the regulation of circadian rhythm. Three transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Feb 2014],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using E4BP4 Monoclonal Antibody