



Eme1 Monoclonal Antibody

Catalog No	YP-mAb-01694
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EME1
Immunogen	Synthesized peptide derived from Eme1 . at AA range: 250-330
Specificity	Eme1 Monoclonal Antibody detects endogenous levels of Eme1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EME1;hMMS4;MMS4;MMS4 homolog;MMS4L;Crossover junction endonuclease EME1
Observed Band	65kD
Calculated Molecular Weight	65kD
Cell Pathway	Nucleus, nucleolus . Recruited to regions of DNA damage in S-phase cells.
Tissue Specificity	Epithelium,Placenta,
Function	cofactor:Magnesium.,function:Interacts with MUS81 to form a DNA structure-specific endonuclease with substrate preference for branched DNA structures with a 5'-end at the branch nick. Typical substrates include 3'-flap structures, replication forks and nicked Holliday junctions. May be required in mitosis for the processing of stalled or collapsed replication forks.,similarity:Belongs to the EME1/MMS4 family.,subcellular location:Recruited to regions of DNA damage in S-phase cells.,subunit:May self-associate. Interacts with MUS81. Interacts with ERCC4 and FANCM.,
Background	This gene encodes a protein that complexes with methyl methanesulfonate-sensitive UV-sensitive 81 protein to form an endonuclease complex. The encoded protein interacts with specific DNA structures including nicked Holliday



junctions, 3'-flap structures and aberrant replication fork structures. This protein may be involved in repairing DNA damage and in maintaining genomic stability. Alternative splicing results in multiple transcript variants.[provided by RefSeq, Oct 2009],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Eme1 Monoclonal Antibody