



ERCC4 Monoclonal Antibody

Catalog No	YP-mAb-01700
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ERCC4
Immunogen	The antiserum was produced against synthesized peptide derived from human XPF. AA range:801-850
Specificity	ERCC4 Monoclonal Antibody detects endogenous levels of ERCC4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ERCC4;DNA repair endonuclease XPF;ERCC11;XPF;DNA excision repair protein ERCC-4;DNA repair protein complementing XP-F cells;Xeroderma pigmentosum group F-complementing protein
Observed Band	103kD
Calculated Molecular Weight	103kD
Cell Pathway	Nucleus . Chromosome . Localizes to sites of DNA damage. .
Tissue Specificity	Epithelium,Fibroblast,
Function	cofactor:Magnesium.,disease:Defects in ERCC4 are a cause of XFE progeroid syndrome [MIM:610965]. This syndrome is illustrated by one patient who presented with dwarfism, cachexia and microcephaly.,disease:Defects in ERCC4 are the cause of xeroderma pigmentosum complementation group F (XP-F) [MIM:278760]; also known as xeroderma pigmentosum VI (XP6). XP-F is an autosomal recessive disease characterized by hypersensitivity of the skin to sunlight followed by high incidence of skin cancer and frequent neurologic abnormalities.,function:Structure-specific DNA repair endonuclease responsible for the 5-prime incision during DNA repair. Involved in homologous recombination that assists in removing interstrand cross-link.,similarity:Belongs



to the XPF family.,subunit:Heterodimer composed of ERCC1 and XPF/ERCC4. Interacts with EME1.,

Background

The protein encoded by this gene forms a complex with ERCC1 and is involved in the 5' incision made during nucleotide excision repair. This complex is a structure specific DNA repair endonuclease that interacts with EME1. Defects in this gene are a cause of xeroderma pigmentosum complementation group F (XP-F), or xeroderma pigmentosum VI (XP6).[provided by RefSeq, Mar 2009],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ERCC4 Monoclonal Antibody