



FAST-1/2 Monoclonal Antibody

Catalog No	YP-mAb-01706
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FOXH1
Immunogen	The antiserum was produced against synthesized peptide derived from human FOXH1. AA range:10-59
Specificity	FAST-1/2 Monoclonal Antibody detects endogenous levels of FAST-1/2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FAST-1/2;FOXH1;Forkhead box protein H1;FAST1;FAST2;Forkhead activin signal transducer 1;Fast-1;hFAST-1;Forkhead activin signal transducer 2;Fast-2
Observed Band	40kD
Calculated Molecular Weight	40kD
Cell Pathway	Nucleus.
Tissue Specificity	Ubiquitous.
Function	domain:The FM region is required for binding smad2/smad4 complexes. FM2 is more effective than FM1 and only interacts with phosphorylated smad2 that is in an activated smad complex.,function:Transcriptional activator. Recognizes and binds to the DNA sequence 5'-TGT[GT][GT]ATT-3'. Required for induction of the goosecoid (GSC) promoter by TGF-beta or activin signaling. Forms a transcriptionally active complex containing FOXH1/SMAD2/SMAD4 on a site on the GSC promoter called TARE (TGF-beta/activin response element).,similarity:Contains 1 fork-head DNA-binding domain.,subunit:Interacts with the MH2 domains of SMAD2 and SMAD3.,tissue specificity:Ubiquitous.,



Background

FOXH1 encodes a human homolog of Xenopus forkhead activin signal transducer-1. FOXH1 protein binds SMAD2 and activates an activin response element via binding the DNA motif TGT(G/T)(T/G)ATT. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using FAST-1/2 Monoclonal Antibody