



FBP3 Monoclonal Antibody

Catalog No	YP-mAb-01708
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FUBP3
Immunogen	The antiserum was produced against synthesized peptide derived from human FUBP3. AA range:201-250
Specificity	FBP3 Monoclonal Antibody detects endogenous levels of FBP3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FBP3;FUBP3;FUSE binding protein 3;Far upstream element-binding protein 3; FUSE-binding protein 3
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Nucleus .
Tissue Specificity	Detected in a number of cell lines.
Function	function:May interact with single-stranded DNA from the far-upstream element (FUSE). May activate gene expression.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 4 KH domains.,tissue specificity:Detected in a number of cell lines.,
Background	function:May interact with single-stranded DNA from the far-upstream element (FUSE). May activate gene expression.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 4 KH domains.,tissue specificity:Detected in a number of cell lines.,



matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using FBP3 Monoclonal Antibody