




FoxB1/2 Monoclonal Antibody

Catalog No	YP-mAb-01717
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FOXB1/FOXB2
Immunogen	The antiserum was produced against synthesized peptide derived from human FOXB1/2. AA range:51-100
Specificity	FoxB1/2 Monoclonal Antibody detects endogenous levels of FoxB1/2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FoxB1/2;FOXB1/FOXB2;Forkhead box protein B1/2;FOXB1;FKH5;Forkhead box protein B1;Transcription factor FKH-5;FOXB2;Forkhead box protein B2
Observed Band	45kD
Calculated Molecular Weight	45kD
Cell Pathway	Nucleus .
Tissue Specificity	Lung,Placenta,
Function	similarity:Contains 1 fork-head DNA-binding domain.,
Background	similarity:Contains 1 fork-head DNA-binding domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using FoxB1/2 Monoclonal Antibody