




FoxB2 Monoclonal Antibody

Catalog No	YP-mAb-01718
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FOXB2
Immunogen	The antiserum was produced against synthesized peptide derived from human FOXB2. AA range:251-300
Specificity	FoxB2 Monoclonal Antibody detects endogenous levels of FoxB2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FoxB2;Forkhead box protein B2
Observed Band	46kD
Calculated Molecular Weight	46kD
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:Transcription factor.,similarity:Contains 1 fork-head DNA-binding domain.,
Background	function:Transcription factor.,similarity:Contains 1 fork-head DNA-binding domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using FoxB2
Monoclonal Antibody