




FoxL1 Monoclonal Antibody

Catalog No	YP-mAb-01729
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FOXL1
Immunogen	The antiserum was produced against synthesized peptide derived from human FOXL1. AA range:121-170
Specificity	FoxL1 Monoclonal Antibody detects endogenous levels of FoxL1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FoxL1;Forkhead box protein L1;FKHL11;FREAC7;Forkhead-related protein FKHL11;Forkhead-related transcription factor 7;FREAC-7
Observed Band	36kD
Calculated Molecular Weight	36kD
Cell Pathway	Nucleus.
Tissue Specificity	Chondrosarcoma Lung Metastasis,
Function	similarity:Contains 1 fork-head DNA-binding domain.,
Background	This gene encodes a member of the forkhead/winged helix-box (FOX) family of transcription factors. FOX transcription factors are characterized by a distinct DNA-binding forkhead domain and play critical roles in the regulation of multiple processes including metabolism, cell proliferation and gene expression during ontogenesis. [provided by RefSeq, Nov 2012],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using FoxL1 Monoclonal Antibody