



FoxN2 Monoclonal Antibody


Catalog No	YP-mAb-01730
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FOXN2
Immunogen	The antiserum was produced against synthesized peptide derived from human FOXN2. AA range:341-390
Specificity	FoxN2 Monoclonal Antibody detects endogenous levels of FoxN2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FoxN2;Forkhead box protein N2;HTLF;Human T-cell leukemia virus enhancer factor
Observed Band	36kD
Calculated Molecular Weight	36kD
Cell Pathway	Nucleus.
Tissue Specificity	PCR rescued clones,Placenta,
Function	function:Binds to the purine-rich region in HTLV-I LTR.,similarity:Contains 1 fork-head DNA-binding domain.,
Background	This gene encodes a forkhead domain binding protein and may function in the transcriptional regulation of the human T-cell leukemia virus long terminal repeat. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using FoxN2 Monoclonal Antibody