



FoxR1 Monoclonal Antibody

Catalog No	YP-mAb-01740
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FOXR1
Immunogen	The antiserum was produced against synthesized peptide derived from human FOXR1. AA range:231-280
Specificity	FoxR1 Monoclonal Antibody detects endogenous levels of FoxR1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FoxR1;Forkhead box protein R1;FOXN5;DLNB13;Forkhead box protein N5
Observed Band	35kD
Calculated Molecular Weight	35kD
Cell Pathway	Cytoplasm . Cytoplasm, perinuclear region . Located in the cytoplasm of spermatocytes and strongly accumulates at the perinuclear region in elongated spermatids. .
Tissue Specificity	Expressed in testis (at protein level).
Function	similarity:Contains 1 fork-head DNA-binding domain.,
Background	This gene encodes a member of the forkhead box (FOX) family of transcription factors. FOX family members are monomeric, helix-turn-helix proteins with a core DNA-binding domain of approximately 110 aa. Many FOX transcription factors play roles in determining cell fates during early development. This forkhead box protein lacks the C-terminal basic region found in many other FOX family members. It is located within the 11q23.3 region which is commonly deleted in neuroblastomas. [provided by RefSeq, Jul 2008],




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using FoxR1 Monoclonal Antibody