



Helios Monoclonal Antibody

Catalog No	YP-mAb-01770
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	IKZF2
Immunogen	The antiserum was produced against synthesized peptide derived from human IKZF2. AA range:331-380
Specificity	Helios Monoclonal Antibody detects endogenous levels of Helios protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IKZF2;ZNFN1A2;ZNF1A2;Zinc finger protein Helios;IKZF 2;Helios;Ikaros family zinc finger protein 2
Observed Band	56kD
Calculated Molecular Weight	56kD
Cell Pathway	Nucleus .
Tissue Specificity	Epithelium,Testis,
Function	function:Associates with Ikaros at centromeric heterochromatin., similarity:Belongs to the Ikaros C2H2-type zinc-finger protein family., similarity:Contains 6 C2H2-type zinc fingers.,subunit:Interacts with IKZF4 AND IKZF5.,
Background	This gene encodes a member of the Ikaros family of zinc-finger proteins. Three members of this protein family (Ikaros, Aiolos and Helios) are hematopoietic-specific transcription factors involved in the regulation of lymphocyte development. This protein forms homo- or hetero-dimers with other Ikaros family members, and is thought to function predominantly in early hematopoietic development. Multiple transcript variants encoding different isoforms have been



found for this gene, but the biological validity of some variants has not been determined. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Helios Monoclonal Antibody