



hnRNP LL Monoclonal Antibody

Catalog No	YP-mAb-01810
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HNRPLL
Immunogen	The antiserum was produced against synthesized peptide derived from human HNRPLL. AA range:241-290
Specificity	hnRNP LL Monoclonal Antibody detects endogenous levels of hnRNP LL protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Stromal RNA-regulating factor;SRRF;hnRNPLL;Heterogeneous nuclear ribonucleoprotein L-like;hnRNP LL;HNRPLL;BLOCK24
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	nucleus,membrane,viral nucleocapsid,intracellular ribonucleoprotein complex,
Tissue Specificity	Widely expressed. Detected in bone marrow stroma cells, skeletal muscle, heart, placenta, pancreas, kidney and lung.
Function	function:May bind RNA and play a role in mRNA processing.,similarity:Contains 3 RRM (RNA recognition motif) domains.,tissue specificity:Widely expressed. Detected in bone marrow stroma cells, skeletal muscle, heart, placenta, pancreas, kidney and lung.,
Background	heterogeneous nuclear ribonucleoprotein L like(HNRNPLL) Homo sapiens HNRNPLL is a master regulator of activation-induced alternative splicing in T cells. In particular, it alters splicing of CD45 (PTPRC; MIM 151460), a tyrosine phosphatase essential for T-cell development and activation (Oberdoerffer et al., 2008 [PubMed 18669861]).[supplied by OMIM, Aug 2008],




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using hnRNP LL Monoclonal Antibody