




HSFX1 Monoclonal Antibody

Catalog No	YP-mAb-01819
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HSFX1
Immunogen	The antiserum was produced against synthesized peptide derived from human LW-1. AA range:10-59
Specificity	HSFX1 Monoclonal Antibody detects endogenous levels of HSFX1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HSFX1;Heat shock transcription factor X-linked;LW-1;HSFX2;Heat shock transcription factor;X-linked
Observed Band	47kD
Calculated Molecular Weight	47kD
Cell Pathway	Nucleus . Cytoplasm .
Tissue Specificity	Testis-specific.
Function	similarity:Belongs to the HSF family.,tissue specificity:Testis-specific.,
Background	similarity:Belongs to the HSF family.,tissue specificity:Testis-specific.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using HSF1
Monoclonal Antibody