



KOX11 Monoclonal Antibody


Catalog No	YP-mAb-01842
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF18
Immunogen	The antiserum was produced against synthesized peptide derived from human ZNF18. AA range:261-310
Specificity	KOX11 Monoclonal Antibody detects endogenous levels of KOX11 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HDSG1;KOX11;Zfp535;zinc finger protein 18;Zinc finger protein KOX11; ZKSCAN6;ZNF18;ZNF535;Heart development-specific gene 1 protein;Zinc finger protein 535;Zinc finger protein with KRAB and SCAN domains 6
Observed Band	62kD
Calculated Molecular Weight	62kD
Cell Pathway	Nucleus .
Tissue Specificity	Brain,Cerebellum,Lymphoid,Prostate,
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 1 SCAN box domain.,similarity:Contains 5 C2H2-type zinc fingers.,
Background	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 1 SCAN box domain.,similarity:Contains 5 C2H2-type zinc fingers.,

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using KOX11 Monoclonal Antibody