



MRRF Monoclonal Antibody

Catalog No	YP-mAb-01877
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MRRF
Immunogen	The antiserum was produced against synthesized peptide derived from human MRRF. AA range:181-230
Specificity	MRRF Monoclonal Antibody detects endogenous levels of MRRF protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RRF;Ribosome-releasing factor;mitochondrial;Ribosome-recycling factor;mtRRF;MRFF;MRRF;Ribosome-recycling factor mitochondrial; mitochondrial
Observed Band	28kD
Calculated Molecular Weight	28kD
Cell Pathway	Mitochondrion .
Tissue Specificity	Lung,Placenta,
Function	function:Responsible for the release of ribosomes from messenger RNA at the termination of protein biosynthesis. May increase the efficiency of translation by recycling ribosomes from one round of translation to another.,similarity:Belongs to the RRF family.,
Background	This gene encodes a ribosome recycling factor, which is a component of the mitochondrial translational machinery. The encoded protein, along with mitochondrial elongation factor 2, functions in ribosomal recycling at the termination of mitochondrial translation by mediating the disassembly of ribosomes from messenger RNA. A pseudogene of this gene has been identified on chromosome X. [provided by RefSeq, Oct 2016],




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using MRRF Monoclonal Antibody