



Myf-6 Monoclonal Antibody

Catalog No	YP-mAb-01886
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MYF6
Immunogen	The antiserum was produced against synthesized peptide derived from human MYF6. AA range:116-165
Specificity	Myf-6 Monoclonal Antibody detects endogenous levels of Myf-6 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Myf-6;Myf 6;Muscle-specific regulatory factor 4;MRF4;Class C basic helix-loop-helix protein 4;MYF6;Myogenic factor 6;BHLHC4
Observed Band	26kD
Calculated Molecular Weight	26kD
Cell Pathway	Nucleus.
Tissue Specificity	Skeletal muscle.
Function	disease:Defects in MYF6 may be a cause of centronuclear myopathy autosomal dominant (ADCNM) [MIM:160150]; also known as autosomal dominant myotubular myopathy. Centronuclear myopathies are congenital muscle disorders characterized by progressive muscular weakness and wasting involving mainly limb girdle, trunk, and neck muscles. It may also affect distal muscles. Weakness may be present during childhood or adolescence or may not become evident until the third decade of life. Ptosis is a frequent clinical feature. The most prominent histopathologic features include high frequency of centrally located nuclei in muscle fibers not secondary to regeneration, radial arrangement of sarcoplasmic strands around the central nuclei, and predominance and hypotrophy of type 1 fibers.,function:Involved in muscle



differentiation (myogenic factor). Induces fibroblasts to differentiate into myoblasts. Pr

Background

myogenic factor 6(MYF6) Homo sapiens The protein encoded by this gene is a probable basic helix-loop-helix (bHLH) DNA binding protein involved in muscle differentiation. The encoded protein likely acts as a heterodimer with another bHLH protein. Defects in this gene are a cause of autosomal dominant centronuclear myopathy (ADCNM). [provided by RefSeq, May 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Myf-6 Monoclonal Antibody