



MZF-1 Monoclonal Antibody

Catalog No	YP-mAb-01888
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MZF1
Immunogen	The antiserum was produced against synthesized peptide derived from human MZF-1. AA range:1-50
Specificity	MZF-1 Monoclonal Antibody detects endogenous levels of MZF-1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MZF-1;MZF1;Myeloid zinc finger 1;MZF;ZNF42;ZSCAN6;Zinc finger and SCAN domain-containing protein 6;Zinc finger protein 42
Observed Band	82kD
Calculated Molecular Weight	82kD
Cell Pathway	Nucleus.
Tissue Specificity	Preferentially expressed in differentiating myeloid cells. Detected in osteoblasts.
Function	function:May be one regulator of transcriptional events during hemopoietic development.,induction:By retinoic acid.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 SCAN box domain.,similarity:Contains 13 C2H2-type zinc fingers.,tissue specificity:Preferentially expressed in differentiating myeloid cells.,
Background	function:May be one regulator of transcriptional events during hemopoietic development.,induction:By retinoic acid.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 SCAN box domain.,similarity:Contains 13 C2H2-type zinc fingers.,tissue specificity:Preferentially expressed in differentiating myeloid cells.,


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using MZF-1 Monoclonal Antibody