



NF-YA Monoclonal Antibody

Catalog No	YP-mAb-01901
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NFYA
Immunogen	The antiserum was produced against synthesized peptide derived from human NFYA. AA range:261-310
Specificity	NF-YA Monoclonal Antibody detects endogenous levels of NF-YA protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NF-YA;NF YA;HAP2;CBF B;CBF A;NFYA;Nuclear transcription factor Y subunit alpha;CAAT box DNA-binding protein subunit A;Nuclear transcription factor Y subunit A
Observed Band	43kD
Calculated Molecular Weight	43kD
Cell Pathway	Nucleus.
Tissue Specificity	Testis,
Function	function:Stimulates the transcription of various genes by recognizing and binding to a CCAAT motif in promoters, for example in type 1 collagen, albumin and beta-actin genes.,similarity:Belongs to the NFYA/HAP2 subunit family., subunit:Heterotrimeric transcription factor composed of three components, NF-YA, NF-YB and NF-YC. NF-YB and NF-YC must interact and dimerize for NF-YA association and DNA binding.,
Background	The protein encoded by this gene is one subunit of a trimeric complex, forming a highly conserved transcription factor that binds to CCAAT motifs in the promoter regions in a variety of genes. Subunit A associates with a tight dimer composed



of the B and C subunits, resulting in a trimer that binds to DNA with high specificity and affinity. The sequence specific interactions of the complex are made by the A subunit, suggesting a role as the regulatory subunit. In addition, there is evidence of post-transcriptional regulation in this gene product, either by protein degradation or control of translation. Further regulation is represented by alternative splicing in the glutamine-rich activation domain, with clear tissue-specific preferences for the two isoforms. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using NF-YA Monoclonal Antibody