



NTPase Monoclonal Antibody


Catalog No	YP-mAb-01921
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NTPCR
Immunogen	The antiserum was produced against synthesized peptide derived from human C1orf57. AA range:141-190
Specificity	NTPase Monoclonal Antibody detects endogenous levels of NTPase protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C1orf57;NTPase;RP4 678E16.2;NTPCR;Cancer-related nucleoside-triphosphatase;Nucleoside triphosphate phosphohydrolase
Observed Band	21kD
Calculated Molecular Weight	21kD
Cell Pathway	membrane,extracellular exosome,
Tissue Specificity	Fetal brain,Skin,
Function	catalytic activity:NTP + H(2)O = NDP + phosphate.,function:Has nucleotide phosphatase activity towards ATP, GTP, CTP, TTP and UTP. Hydrolyzes nucleoside diphosphates with lower efficiency.,similarity:Belongs to the THEP1 NTPase family.,subunit:Monomer.,
Background	The protein encoded by this gene is a non-specific nucleoside triphosphatase that is slow-acting in vitro. This gene is overexpressed in many tumor tissues, and while it is not essential for the cell, overexpression is cytotoxic. However, the cytotoxicity is not related to its triphosphatase activity. [provided by RefSeq, Jul 2016],

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using NTPase Monoclonal Antibody