



Per3 Monoclonal Antibody

Catalog No	YP-mAb-01938
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PER3
Immunogen	The antiserum was produced against synthesized peptide derived from human PER3. AA range:21-70
Specificity	Per3 Monoclonal Antibody detects endogenous levels of Per3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Cell growth-inhibiting gene 13 protein;GIG13;hPER3;period homolog 3 (Drosophila);Per3;Period circadian protein homolog 3;Circadian clock protein PERIOD 3
Observed Band	43kD
Calculated Molecular Weight	132KD
Cell Pathway	Cytoplasm . Nucleus . Mainly cytoplasmic. Translocates to the nucleus through binding PER1, PER2, CRY1 or CRY2, but not TIMELESS. .
Tissue Specificity	Brain,
Function	function:Component of the circadian clock mechanism which is essential for generating circadian rhythms. Function unknown.,PTM:Phosphorylation appears to require association with PER1 and translocation to the nucleus.,PTM:Ubiquitinated.,similarity:Contains 1 PAC (PAS-associated C-terminal) domain.,similarity:Contains 2 PAS (PER-ARNT-SIM) domains.,subcellular location:Mainly cytoplasmic. Translocates to the nucleus through binding PER1, PER2, CRY1 or CRY2, but not TIMELESS.,subunit:Homodimer. Component of the circadian core oscillator, which includes the CRY proteins, CLOCK or NPAS2, BMAL1 or BMAL2, CSNK1D and/or CSNK1E, TIMELESS and the PER



proteins. Interacts directly with PER1, PER2, CRY1, CRY2, and TIMELESS. Interaction with CSNK1D or CSNK1E promotes nuclear location of PER proteins.,

Background

This gene is a member of the Period family of genes and is expressed in a circadian pattern in the suprachiasmatic nucleus, the primary circadian pacemaker in the mammalian brain. Genes in this family encode components of the circadian rhythms of locomotor activity, metabolism, and behavior. This gene is upregulated by CLOCK/ARNTL heterodimers but then represses this upregulation in a feedback loop using PER/CRY heterodimers to interact with CLOCK/ARNTL. Polymorphisms in this gene have been linked to sleep disorders. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2014],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Per3 Monoclonal Antibody