



RBAK Monoclonal Antibody

Catalog No	YP-mAb-01967
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RBAK
Immunogen	The antiserum was produced against synthesized peptide derived from human RBAK. AA range:1-50
Specificity	RBAK Monoclonal Antibody detects endogenous levels of RBAK protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RBAK;RB-associated KRAB zinc finger protein;ZNF769;RB-associated KRAB repressor;hRBaK;Zinc finger protein 769
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in bone, brain, heart, kidney, liver, lung, pancreas and placenta.
Function	function:May repress E2F-dependent transcription. May promote AR-dependent transcription.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 16 C2H2-type zinc fingers.,subunit:Interacts with AR and RB1. May also interact with other nuclear hormone receptors such as NR3C1/GR.,tissue specificity:Expressed in bone, brain, heart, kidney, liver, lung, pancreas and placenta.,
Background	This gene encodes a nuclear protein which interacts with the tumor suppressor retinoblastoma 1. The two interacting proteins are thought to act as a transcriptional repressor for promoters which are activated by the E2F1 transcription factor. This protein contains a Kruppel-associated box (KRAB),



which is a transcriptional repressor motif. Read-through transcripts that include exons from the downstream gene LOC389458 are expressed from this locus. [provided by RefSeq, Mar 2011],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using RBAK Monoclonal Antibody