




RBM16 Monoclonal Antibody

Catalog No	YP-mAb-01968
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SCAF8
Immunogen	Synthesized peptide derived from RBM16 . at AA range: 550-630
Specificity	RBM16 Monoclonal Antibody detects endogenous levels of RBM16 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SCAF8;RNA-binding motif protein 16;RNA binding protein 16;KIAA1116;CDC5L complex-associated protein 7;RBM16;Protein SCAF8;CCAP7;SR-related and CTD-associated factor 8
Observed Band	140kD
Calculated Molecular Weight	140kD
Cell Pathway	Nucleus . Nucleus matrix . Detected in granular nuclear foci which correspond to sites of active transcription. .
Tissue Specificity	Brain,Epithelium,Testis,Thyroid gland,Uterus,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 CID domain.,similarity:Contains 1 RRM (RNA recognition motif) domain.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 CID domain.,similarity:Contains 1 RRM (RNA recognition motif) domain.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using RBM16 Monoclonal Antibody