




RBM34 Monoclonal Antibody

Catalog No	YP-mAb-01970
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RBM34
Immunogen	The antiserum was produced against synthesized peptide derived from human RBM34. AA range:345-394
Specificity	RBM34 Monoclonal Antibody detects endogenous levels of RBM34 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA0117;RNA-binding motif protein 34;RNA-binding protein 34;RBM34
Observed Band	50kD
Calculated Molecular Weight	50kD
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	Bone marrow,Epithelium,Lung,Testis,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Belongs to the RRM RBM34 family.,similarity:Contains 2 RRM (RNA recognition motif) domains.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Belongs to the RRM RBM34 family.,similarity:Contains 2 RRM (RNA recognition motif) domains.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using RBM34 Monoclonal Antibody