




RECS1 Monoclonal Antibody

Catalog No	YP-mAb-01972
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMBIM1
Immunogen	The antiserum was produced against synthesized peptide derived from human TMBIM1. AA range:251-300
Specificity	RECS1 Monoclonal Antibody detects endogenous levels of RECS1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RECS1;TMBIM1;Protein lifeguard 3;LFG3;PP1201;PSEC0158;Protein RECS1 homolog;Transmembrane BAX inhibitor motif-containing protein 1
Observed Band	35kD
Calculated Molecular Weight	35kD
Cell Pathway	Membrane ; Multi-pass membrane protein . Lysosome membrane . Endosome membrane .
Tissue Specificity	Aorta,Brain,Colon,Pancreas,Placenta,Platelet,Tongue,Trachea
Function	similarity:Belongs to the BI1 family.,
Background	similarity:Belongs to the BI1 family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using RECS1 Monoclonal Antibody