



RFWD3 Monoclonal Antibody

Catalog No	YP-mAb-01978
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RFWD3
Immunogen	The antiserum was produced against synthesized peptide derived from human RFWD3. AA range:374-423
Specificity	RFWD3 Monoclonal Antibody detects endogenous levels of RFWD3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RFWD3;E3 ubiquitin-protein ligase RFWD3;RNF201;RING finger and WD repeat domain-containing protein 3;RING finger protein 201
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Nucleus . Nucleus, PML body . Cytoplasm . In undamaged cells, found both in the cytoplasm and in the nucleus, partially associated with PML nuclear bodies (PubMed:21558276). In response to replication block, such as that caused by hydroxyurea treatment, or to DNA damage caused by ionizing radiations or doxorubicin, recruited to the nucleus, to stalled replication forks or to sites of DNA repair (PubMed:21504906, PubMed:28575657). This recruitment depends upon RPA2 (PubMed:21504906). .
Tissue Specificity	Brain, Teratocarcinoma, Testis,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 RING-type zinc finger., similarity:Contains 3 WD repeats.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 RING-type zinc finger., similarity:Contains 3 WD repeats.,


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using RFWD3 Monoclonal Antibody