




RGAG1 Monoclonal Antibody

Catalog No	YP-mAb-01981
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RGAG1
Immunogen	The antiserum was produced against synthesized peptide derived from human RGAG1. AA range:1331-1380
Specificity	RGAG1 Monoclonal Antibody detects endogenous levels of RGAG1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RGAG1;Retrotransposon gag domain-containing protein 1;KIAA1318;Tumor antigen BJ-HCC-23
Observed Band	144kD
Calculated Molecular Weight	144kD
Cell Pathway	
Tissue Specificity	Brain,Carcinoma,
Function	
Background	RGAG1 (Retrotransposon Gag Domain Containing 1) is a Protein Coding gene. Diseases associated with RGAG1 include non-syndromic x-linked intellectual disability.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using RGAG1 Monoclonal Antibody