



Ribosomal Protein L10L Monoclonal Antibody

Catalog No	YP-mAb-01982
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RPL10L
Immunogen	The antiserum was produced against synthesized peptide derived from human RPL10L. AA range:151-200
Specificity	Ribosomal Protein L10L Monoclonal Antibody detects endogenous levels of Ribosomal Protein L10L protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	60S ribosomal protein L10 like;ribosomal protein L10 like;RPL10_5_1358; RPL10L;Ribosomal Protein L10L;60S ribosomal protein L10-like
Observed Band	20kD
Calculated Molecular Weight	20kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Almost testis-specific (PubMed:12490704, PubMed:19333399, PubMed:28502657). Also expressed in pre- and postmenopausal ovary (PubMed:19333399).
Function	function:May play a role in compensating for the inactivated X-linked gene during spermatogenesis.,miscellaneous:This gene has no introns in its coding regions, and therefore was most likely produced by retrotransposition of the original X-linked gene during evolution.,similarity:Belongs to the ribosomal protein L10e family.,tissue specificity:Testis specific.,
Background	This gene encodes a protein sharing sequence similarity with ribosomal protein L10. It is not currently known whether the encoded protein is a functional ribosomal protein or whether it has evolved a function that is independent of the



ribosome. This gene is intronless. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Ribosomal Protein L10L Monoclonal Antibody