




RNF113B Monoclonal Antibody

Catalog No	YP-mAb-01985
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RNF113B
Immunogen	Synthesized peptide derived from the Internal region of human RNF113B.
Specificity	RNF113B Monoclonal Antibody detects endogenous levels of RNF113B protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RNF113B;RING finger protein 113B;RNF161;ZNF183L1;Zinc finger protein 183-like 1
Observed Band	36kD
Calculated Molecular Weight	36kD
Cell Pathway	U2-type spliceosomal complex,
Tissue Specificity	Testis,
Function	similarity:Contains 1 C3H1-type zinc finger.,similarity:Contains 1 RING-type zinc finger.,
Background	similarity:Contains 1 C3H1-type zinc finger.,similarity:Contains 1 RING-type zinc finger.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using RNF113B Monoclonal Antibody