



RNF15 Monoclonal Antibody

Catalog No	YP-mAb-01989
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TRIM38
Immunogen	The antiserum was produced against synthesized peptide derived from human TRIM38. AA range:161-210
Specificity	RNF15 Monoclonal Antibody detects endogenous levels of RNF15 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TRIM;RORET;RNF15;RING finger protein 15;EC:2.3.2.27;TRIM38;Tripartite motif-containing protein 38;Zinc finger protein RoRet
Observed Band	58kD
Calculated Molecular Weight	58kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Ubiquitous.
Function	similarity:Contains 1 B box-type zinc finger.,similarity:Contains 1 B30.2/SPRY domain.,similarity:Contains 1 RING-type zinc finger.,tissue specificity:Ubiquitous.,
Background	This gene encodes a member of the tripartite motif (TRIM) family. The encoded protein contains a RING-type zinc finger, B box-type zinc finger and SPRY domain. The function of this protein has not been identified. A pseudogene of this gene is located on the long arm of chromosome 4. [provided by RefSeq, Jul 2012],


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using RNF15 Monoclonal Antibody