




SAF-A2 Monoclonal Antibody

Catalog No	YP-mAb-02001
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HNRNPUL2
Immunogen	The antiserum was produced against synthesized peptide derived from human HNRNPUL2. AA range:601-650
Specificity	SAF-A2 Monoclonal Antibody detects endogenous levels of SAF-A2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HNRNPUL2;Heterogeneous nuclear ribonucleoprotein U-like protein 2;HNRNP;HNRNPUL2;SAF A2;SAF-A2;Scaffold-attachment factor A2
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Nucleus .
Tissue Specificity	Epithelium,Liver,Melanoma,Ovarian carcinoma,
Function	similarity:Contains 1 B30.2/SPRY domain.,similarity:Contains 1 SAP domain.,subunit:Binds to MLF1 and retains it in the nucleus.,
Background	similarity:Contains 1 B30.2/SPRY domain.,similarity:Contains 1 SAP domain.,subunit:Binds to MLF1 and retains it in the nucleus.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using SAF-A2 Monoclonal Antibody