



SAP 49 Monoclonal Antibody

Catalog No	YP-mAb-02003
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SF3B4
Immunogen	The antiserum was produced against synthesized peptide derived from human SF3B4. AA range:61-110
Specificity	SAP 49 Monoclonal Antibody detects endogenous levels of SAP 49 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SAP 49;SAP49;SF3B4;SF3b49;SF3b50;Splicing factor 3B subunit 4;Pre-mRNA-splicing factor SF3b 49 kDa subunit;Spliceosome-associated protein 49
Observed Band	55kD
Calculated Molecular Weight	55kD
Cell Pathway	Nucleus .
Tissue Specificity	Ovarian carcinoma,Skin,Testis,Thymus,
Function	function:Subunit of the splicing factor SF3B required for 'A' complex assembly formed by the stable binding of U2 snRNP to the branchpoint sequence (BPS) in pre-mRNA. Sequence independent binding of SF3A/SF3B complex upstream of the branch site is essential, it may anchor U2 snRNP to the pre-mRNA. May also be involved in the assembly of the 'E' complex. SF3B4 has been found in complex 'B' and 'C' as well. Belongs also to the minor U12-dependent spliceosome, which is involved in the splicing of rare class of nuclear pre-mRNA intron.,similarity:Belongs to the SF3B4 family.,similarity:Contains 2 RRM (RNA recognition motif) domains.,subunit:Component of splicing factor SF3B which is composed of at least eight subunits; SF3B1/SAP155/SF3B155, SF3B2/SAP145/SF3B145, SF3B3/SAP130/SF3B130, SF3B4/SAP49/SF3B49, SF3B14A,



PHF5A/SF3B14B, SF3B10 and SF3B125. SF3B associates with the splicing factor SF3A

Background

This gene encodes one of four subunits of the splicing factor 3B. The protein encoded by this gene cross-links to a region in the pre-mRNA immediately upstream of the branchpoint sequence in pre-mRNA in the prespliceosomal complex A. It also may be involved in the assembly of the B, C and E spliceosomal complexes. In addition to RNA-binding activity, this protein interacts directly and highly specifically with subunit 2 of the splicing factor 3B. This protein contains two N-terminal RNA-recognition motifs (RRMs), consistent with the observation that it binds directly to pre-mRNA. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using SAP 49 Monoclonal Antibody