



SIP1 Monoclonal Antibody

Catalog No	YP-mAb-02012
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZEB2
Immunogen	The antiserum was produced against synthesized peptide derived from human ZEB2. AA range:71-120
Specificity	SIP1 Monoclonal Antibody detects endogenous levels of SIP1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA0569;SIP 1;SIP1;Smad interacting protein 1;SMADIP1;ZEB2;ZFHX1B; ZFX1B;Zinc finger E-box-binding homeobox 2;HRIHFB2411;Smad-interacting protein 1;Zinc finger homeobox protein 1b
Observed Band	157kD
Calculated Molecular Weight	157kD
Cell Pathway	Nucleus . Chromosome .
Tissue Specificity	Brain,Fetal brain,
Function	disease:Defects in ZEB2 are the cause of Hirschsprung disease-mental retardation syndrome (Hirschsprung disease) [MIM:235730]; also known as Mowat-Wilson syndrome (MWS). Hirschsprung disease is a rare autosomal dominant complex developmental disorder. Individuals with functional null mutations present with mental retardation, delayed motor development, epilepsy, and a wide spectrum of clinically heterogeneous features suggestive of neurocristopathies at the cephalic, cardiac, and vagal levels. Affected patients show an easily recognizable facial appearance with deep set eyes and hypertelorism, medially divergent, broad eyebrows, prominent columella, pointed chin and uplifted, notched ear lobes. Additionally, the phenotypic spectrum of



facultative congenital anomalies includes short stature, microcephaly, Hirschsprung disease, malformations of the brain (agenesis of corpus callosum, cereb

Background

The protein encoded by this gene is a member of the Zfh1 family of 2-handed zinc finger/homeodomain proteins. It is located in the nucleus and functions as a DNA-binding transcriptional repressor that interacts with activated SMADs. Mutations in this gene are associated with Hirschsprung disease/Mowat-Wilson syndrome. Alternatively spliced transcript variants have been found for this gene.[provided by RefSeq, Jan 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using SIP1 Monoclonal Antibody