



Snrp116 Monoclonal Antibody

Catalog No	YP-mAb-02028
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	EFTUD2
Immunogen	The antiserum was produced against synthesized peptide derived from human EFTUD2. AA range:321-370
Specificity	Snrp116 Monoclonal Antibody detects endogenous levels of Snrp116 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Snu114;SNRNP116;KIAA0031;hSNU114;Elongation factor Tu GTP-binding domain-containing protein 2;EFTUD2;116 kDa U5 small nuclear ribonucleoprotein component;SNU114 homolog;U5 snRNP-specific protein;116 kDa;U5-116 kDa
Observed Band	109kD
Calculated Molecular Weight	109kD
Cell Pathway	Nucleus .
Tissue Specificity	Bone marrow,Epithelium,Muscle,Thalamus,Uterus,
Function	function:Component of the U5 snRNP complex required for pre-mRNA splicing., PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Belongs to the GTP-binding elongation factor family. EF-G/EF-2 subfamily.,subunit:Identified in the spliceosome C complex, at least composed of AQR, ASCC3L1, C19orf29, CDC40, CDC5L, CRNKL1, DDX23, DDX41, DDX48, DDX5, DGCR14, DHX35, DHX38, DHX8, EFTUD2, FRG1, GPATC1, HNRPA1, HNRPA2B1, HNRPA3, HNRPC, HNRPF, HNRPH1, HNRPK, HNRPM, HNRPR, HNRPU, KIAA1160, KIAA1604, LSM2, LSM3, MAGOH, MORG1, MABPC1, PLRG1, PNN, PPIE, PPIL1, PPIL3, PPWD1, PRPF19, PRPF4B,



PRPF6, PRPF8, RALY, RBM22, RBM8A, RBMX, SART1, SF3A1, SF3A2, SF3A3, SF3B1, SF3B2, SF3B3, SFRS1, SKIV2L2, SNRPA1, SNRPB, SNRPB2, SNRPD1, SNRPD2, SNRPD3, SNRPE, SNRPF, SNRPG, SNW1, SRRM1, SRRM2, SYF2, SYNCRIP, TFIP11, THOC4, U2AF1, WDR57, XAB2 and ZCCHC8. Binds GTP. Interacts with PRPF8.,

Background

This gene encodes a GTPase which is a component of the spliceosome complex which processes precursor mRNAs to produce mature mRNAs. Mutations in this gene are associated with mandibulofacial dysostosis with microcephaly. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2012],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Snrp116 Monoclonal Antibody