




# SRE-ZBP Monoclonal Antibody

|                                    |                                                                                                                                                                               |
|------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-mAb-02041                                                                                                                                                                  |
| <b>Isotype</b>                     | IgG                                                                                                                                                                           |
| <b>Reactivity</b>                  | Human;Mouse;Rat                                                                                                                                                               |
| <b>Applications</b>                | WB                                                                                                                                                                            |
| <b>Gene Name</b>                   | ZNF187                                                                                                                                                                        |
| <b>Immunogen</b>                   | The antiserum was produced against synthesized peptide derived from human ZNF187. AA range:51-100                                                                             |
| <b>Specificity</b>                 | SRE-ZBP Monoclonal Antibody detects endogenous levels of SRE-ZBP protein.                                                                                                     |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                       |
| <b>Source</b>                      | Monoclonal, Mouse,IgG                                                                                                                                                         |
| <b>Purification</b>                | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.                                                          |
| <b>Dilution</b>                    | WB 1:500-1:2000                                                                                                                                                               |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                                                       |
| <b>Purity</b>                      | ≥90%                                                                                                                                                                          |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                                                                  |
| <b>Synonyms</b>                    | SRE-ZBP;ZNF187;Zinc finger protein 187;ZSCAN26;Protein SRE-ZBP;Zinc finger and SCAN domain-containing protein 26                                                              |
| <b>Observed Band</b>               | 55kD                                                                                                                                                                          |
| <b>Calculated Molecular Weight</b> | 55kD                                                                                                                                                                          |
| <b>Cell Pathway</b>                | Nucleus .                                                                                                                                                                     |
| <b>Tissue Specificity</b>          | Cervix,Kidney,                                                                                                                                                                |
| <b>Function</b>                    | function:May be involved in transcriptional regulation.,induction:By serum stimulation.,similarity:Contains 1 SCAN box domain.,similarity:Contains 8 C2H2-type zinc fingers., |
| <b>Background</b>                  | function:May be involved in transcriptional regulation.,induction:By serum stimulation.,similarity:Contains 1 SCAN box domain.,similarity:Contains 8 C2H2-type zinc fingers., |
| <b>matters needing attention</b>   | Avoid repeated freezing and thawing!                                                                                                                                          |

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

**Products Images**

Western Blot analysis of various cells using SRE-ZBP Monoclonal Antibody