



SRp20 Monoclonal Antibody

Catalog No	YP-mAb-02043
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SRSF3
Immunogen	The antiserum was produced against synthesized peptide derived from human SFRS3. AA range:111-160
Specificity	SRp20 Monoclonal Antibody detects endogenous levels of SRp20 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SFRS3;Pre-mRNA-splicing factor SRP2;Serine/arginine-rich splicing factor 3; SFRS 3;Splicing factor;arginine/serine-rich 3;SRp20;SRSF3;Pre-mRNA-splicing factor SRP20
Observed Band	25kD
Calculated Molecular Weight	25kD
Cell Pathway	Nucleus . Nucleus speckle . Cytoplasm . Recruited to nuclear speckles following interaction with YTHDC1. .
Tissue Specificity	Brain,Epithelium,Placenta,Tongue,
Function	function:May be involved in RNA processing in relation with cellular proliferation and/or maturation.,PTM:Extensively phosphorylated on serine residues in the RS domain.,similarity:Belongs to the splicing factor SR family.,similarity:Contains 1 RRM (RNA recognition motif) domain.,subunit:Interacts with CPSF6, RBMY1A1 and SFRS12.,
Background	The protein encoded by this gene is a member of the serine/arginine (SR)-rich family of pre-mRNA splicing factors, which constitute part of the spliceosome. Each of these factors contains an RNA recognition motif (RRM) for binding RNA



and an RS domain for binding other proteins. The RS domain is rich in serine and arginine residues and facilitates interaction between different SR splicing factors. In addition to being critical for mRNA splicing, the SR proteins have also been shown to be involved in mRNA export from the nucleus and in translation. Two transcript variants, one protein-coding and the other non-coding, have been found for this gene. [provided by RefSeq, Sep 2010],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using SRp20 Monoclonal Antibody