




SUHW1 Monoclonal Antibody

Catalog No	YP-mAb-02059
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF280A
Immunogen	The antiserum was produced against synthesized peptide derived from human ZNF280A. AA range:441-490
Specificity	SUHW1 Monoclonal Antibody detects endogenous levels of SUHW1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Zinc finger protein 636;Zinc finger protein 280A;Zinc Finger Protein 280;Suppressor of hairy wing homolog 1;SUHW1;ZNF280A;ZNF280;ZNF636;3'OY11.1
Observed Band	50kD
Calculated Molecular Weight	50kD
Cell Pathway	Nucleus .
Tissue Specificity	Melanoma,Testis,
Function	function:May function as a transcription factor.,similarity:Contains 4 C2H2-type zinc fingers.,
Background	This gene encodes a zinc finger protein. The encoded protein contains 4 C2H2-type zinc fingers, which are commonly found in transcription factors. A variety of functions may be performed by this type of zinc finger protein, including the binding of DNA or RNA. [provided by RefSeq, Apr 2014],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using SUHW1 Monoclonal Antibody