



TBX1 Monoclonal Antibody

Catalog No	YP-mAb-02076
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TBX1
Immunogen	The antiserum was produced against synthesized peptide derived from human TBX1. AA range:311-360
Specificity	TBX1 Monoclonal Antibody detects endogenous levels of TBX1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TBX1;T-box transcription factor TBX1;T-box protein 1;Testis-specific T-box protein
Observed Band	43kD
Calculated Molecular Weight	43kD
Cell Pathway	Nucleus .
Tissue Specificity	Skeletal muscle,Testis,
Function	disease:Defects in TBX1 are a cause of conotruncal heart malformations (CTHM) [MIM:217095]. CTHM consist of cardiac outflow tract defects, such as tetralogy of Fallot, pulmonary atresia, double-outlet right ventricle, truncus arteriosus communis, and aortic arch anomalies.,disease:Defects in TBX1 are a cause of DiGeorge syndrome (DGS) [MIM:188400].,disease:Defects in TBX1 are a cause of velocardiofacial syndrome (VCFS) [MIM:192430]., disease:Haploinsufficiency of the TBX1 gene is responsible for most of the physical malformations present in DiGeorge syndrome (DGS) and velocardiofacial syndrome (VCFS) [MIM:188400, 192430]. DGS is characterized by the association of several malformations: hypoplastic thymus and parathyroid glands, congenital conotruncal cardiopathy, and a subtle but characteristic facial



dysmorphology. VCFS is marked by the association of congenital conotruncal heart defect

Background

This gene is a member of a phylogenetically conserved family of genes that share a common DNA-binding domain, the T-box. T-box genes encode transcription factors involved in the regulation of developmental processes. This gene product shares 98% amino acid sequence identity with the mouse ortholog. DiGeorge syndrome (DGS)/velocardiofacial syndrome (VCFS), a common congenital disorder characterized by neural-crest-related developmental defects, has been associated with deletions of chromosome 22q11.2, where this gene has been mapped. Studies using mouse models of DiGeorge syndrome suggest a major role for this gene in the molecular etiology of DGS/VCFS. Several alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using TBX1 Monoclonal Antibody