




TCEAL6 Monoclonal Antibody

Catalog No	YP-mAb-02087
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TCEAL6
Immunogen	The antiserum was produced against synthesized peptide derived from human TCEAL6. AA range:61-110
Specificity	TCEAL6 Monoclonal Antibody detects endogenous levels of TCEAL6 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TCEAL6;Transcription elongation factor A protein-like 6;TCEA-like protein 6; Transcription elongation factor S-II protein-like 6
Observed Band	25kD
Calculated Molecular Weight	25kD
Cell Pathway	Nucleus .
Tissue Specificity	Brain,
Function	function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the TFS-II family. TFA subfamily.,
Background	function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the TFS-II family. TFA subfamily.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using TCEAL6 Monoclonal Antibody