




THAP4 Monoclonal Antibody

Catalog No	YP-mAb-02113
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	THAP4
Immunogen	The antiserum was produced against synthesized peptide derived from human THAP4. AA range:341-390
Specificity	THAP4 Monoclonal Antibody detects endogenous levels of THAP4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	THAP4;THAP domain-containing protein 4;CGI-36;PP238
Observed Band	63kD
Calculated Molecular Weight	63kD
Cell Pathway	Cytoplasm. Nucleus. Localizes mainly in the cytoplasm and partially in the nucleus. .
Tissue Specificity	Brain,Eye,Melanoma,Muscle,Ovarian carcinoma,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 THAP-type zinc finger.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Contains 1 THAP-type zinc finger.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using THAP4 Monoclonal Antibody