



ZAC1 Monoclonal Antibody

Catalog No	YP-mAb-02161
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PLAGL1
Immunogen	The antiserum was produced against synthesized peptide derived from human PLAGL1. AA range:311-360
Specificity	ZAC1 Monoclonal Antibody detects endogenous levels of ZAC1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Tumor suppressor ZAC;Pleiomorphic adenoma-like protein 1;LOT-1;LOT1;LOT 1;ZAC1;PLAGL1;Zinc finger protein PLAGL1;ZAC;Lost on transformation 1; Tumor supressor ZAC
Observed Band	50kD
Calculated Molecular Weight	50kD
Cell Pathway	Nucleus .
Tissue Specificity	Colon,Fetal kidney,Lung,Mammary gland,Pituitary,Placenta,
Function	function:Shows weak transcriptional activatory activity. Transcriptional regulator of the type 1 receptor for pituitary adenylate cyclase-activating polypeptide., similarity:Belongs to the krueppel C2H2-type zinc-finger protein family., similarity:Contains 7 C2H2-type zinc fingers.,
Background	This gene encodes a C2H2 zinc finger protein that functions as a suppressor of cell growth. This gene is often deleted or methylated and silenced in cancer cells. In addition, overexpression of this gene during fetal development is thought to be the causal factor for transient neonatal diabetes mellitus (TNDM). Alternative splicing and the use of alternative promoters results in multiple



transcript variants encoding two different protein isoforms. The P1 downstream promoter of this gene is imprinted, with preferential expression from the paternal allele in many tissues. [provided by RefSeq, Nov 2015],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ZAC1 Monoclonal Antibody