



ZBP-89 Monoclonal Antibody

Catalog No	YP-mAb-02162
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF148
Immunogen	The antiserum was produced against synthesized peptide derived from human ZNF148. AA range:61-110
Specificity	ZBP-89 Monoclonal Antibody detects endogenous levels of ZBP-89 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BERF 1;BFCOL1;HT BETA;pHZ 52;Transcription factor ZBP 89;ZBP-89;ZNF148;Zinc finger protein 148;ZBP89;Transcription factor ZBP-89;Zinc finger DNA-binding protein 89
Observed Band	89kD
Calculated Molecular Weight	89kD
Cell Pathway	Nucleus.
Tissue Specificity	Epithelium,Insulinoma,Kidney,Peripheral blood,Primary B-Cells,Testis,
Function	function:Involved in transcriptional regulation. Represses the transcription of a number of genes including gastrin, stromelysin and enolase. Binds to the G-rich box in the enhancer region of these genes.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 4 C2H2-type zinc fingers.,subunit:Interacts with HNRPD.,
Background	function:Involved in transcriptional regulation. Represses the transcription of a number of genes including gastrin, stromelysin and enolase. Binds to the G-rich box in the enhancer region of these genes.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 4 C2H2-type zinc



fingers.,subunit:Interacts with HNRPDL.,


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ZBP-89 Monoclonal Antibody