



ZNF287 Monoclonal Antibody

Catalog No	YP-mAb-02182
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF287
Immunogen	The antiserum was produced against synthesized peptide derived from human ZNF287. AA range:241-290
Specificity	ZNF287 Monoclonal Antibody detects endogenous levels of ZNF287 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	zinc finger protein 287;ZKSCAN13;ZNF287;Zinc finger protein with KRAB and SCAN domains 13
Observed Band	90kD
Calculated Molecular Weight	90kD
Cell Pathway	Nucleus .
Tissue Specificity	Brain,
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 1 SCAN box domain.,similarity:Contains 14 C2H2-type zinc fingers.,
Background	zinc finger protein 287(ZNF287) Homo sapiens This gene encodes a member of the krueppel family of zinc finger proteins, suggesting a role as a transcription factor. Its specific function has not been determined. This gene is located near the Smith-Magenis syndrome region on chromosome 17. [provided by RefSeq, Jul 2008],




matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ZNF287 Monoclonal Antibody