




ZNF337 Monoclonal Antibody

Catalog No	YP-mAb-02185
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF337
Immunogen	The antiserum was produced against synthesized peptide derived from human ZNF337. AA range:471-520
Specificity	ZNF337 Monoclonal Antibody detects endogenous levels of ZNF337 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZNF337;Zinc finger protein 337
Observed Band	86kD
Calculated Molecular Weight	86kD
Cell Pathway	Nucleus .
Tissue Specificity	Brain,Lung,
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 20 C2H2-type zinc fingers.,
Background	zinc finger protein 337(ZNF337) Homo sapiens This gene encodes a zinc finger domain containing protein. The function of this protein has yet to be determined. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2014],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ZNF337 Monoclonal Antibody