



ZNF436 Monoclonal Antibody


Catalog No	YP-mAb-02188
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF436
Immunogen	The antiserum was produced against synthesized peptide derived from human ZNF436. AA range:71-120
Specificity	ZNF436 Monoclonal Antibody detects endogenous levels of ZNF436 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZNF436;Zinc finger protein 436;KIAA1710
Observed Band	50kD
Calculated Molecular Weight	50kD
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in fetal brain, heart, liver, spleen, bladder, lung, skin, skeletal muscle, stomach and pancreas.
Function	function:May be a transcriptional repressor.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 12 C2H2-type zinc fingers.,tissue specificity:Expressed in fetal brain, heart, liver, spleen, bladder, lung, skin, skeletal muscle, stomach and pancreas.,
Background	function:May be a transcriptional repressor.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 12 C2H2-type zinc fingers.,tissue specificity:Expressed in fetal brain, heart, liver, spleen, bladder, lung, skin, skeletal muscle, stomach and pancreas.,

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ZNF436 Monoclonal Antibody