




ZNF596 Monoclonal Antibody

Catalog No	YP-mAb-02195
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF596
Immunogen	The antiserum was produced against synthesized peptide derived from human ZNF596. AA range:261-310
Specificity	ZNF596 Monoclonal Antibody detects endogenous levels of ZNF596 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZNF596;Zinc finger protein 596
Observed Band	58kD
Calculated Molecular Weight	58kD
Cell Pathway	Nucleus .
Tissue Specificity	Testis,
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 11 C2H2-type zinc fingers.,
Background	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 11 C2H2-type zinc fingers.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western Blot analysis of various cells using ZNF596 Monoclonal Antibody