




# ZNF608 Monoclonal Antibody

|                                    |                                                                                                                                                                                                          |
|------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-mAb-02197                                                                                                                                                                                             |
| <b>Isotype</b>                     | IgG                                                                                                                                                                                                      |
| <b>Reactivity</b>                  | Human;Mouse;Rat                                                                                                                                                                                          |
| <b>Applications</b>                | WB                                                                                                                                                                                                       |
| <b>Gene Name</b>                   | ZNF608                                                                                                                                                                                                   |
| <b>Immunogen</b>                   | The antiserum was produced against synthesized peptide derived from human ZNF608. AA range:1231-1280                                                                                                     |
| <b>Specificity</b>                 | ZNF608 Monoclonal Antibody detects endogenous levels of ZNF608 protein.                                                                                                                                  |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                  |
| <b>Source</b>                      | Monoclonal, Mouse,IgG                                                                                                                                                                                    |
| <b>Purification</b>                | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                     |
| <b>Dilution</b>                    | WB 1:500-1:2000                                                                                                                                                                                          |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                                                                                                  |
| <b>Purity</b>                      | ≥90%                                                                                                                                                                                                     |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                                                                                             |
| <b>Synonyms</b>                    | KIAA1281;NY REN 36;zinc finger protein 608;ZNF608;ZNF608-Specific;Renal carcinoma antigen NY-REN-36                                                                                                      |
| <b>Observed Band</b>               | 160kD                                                                                                                                                                                                    |
| <b>Calculated Molecular Weight</b> | 160kD                                                                                                                                                                                                    |
| <b>Cell Pathway</b>                |                                                                                                                                                                                                          |
| <b>Tissue Specificity</b>          | Brain,PNS,Renal cell carcinoma,Skin,Testis,                                                                                                                                                              |
| <b>Function</b>                    | miscellaneous:Antigen recognized by autologous antibody in patients with renal-cell carcinoma.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 C2H2-type zinc finger., |
| <b>Background</b>                  | miscellaneous:Antigen recognized by autologous antibody in patients with renal-cell carcinoma.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 C2H2-type zinc finger., |
| <b>matters needing attention</b>   | Avoid repeated freezing and thawing!                                                                                                                                                                     |

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

**Products Images**

Western Blot analysis of various cells using ZNF608 Monoclonal Antibody