



ZNF785 Monoclonal Antibody

Catalog No	YP-mAb-02205
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF785
Immunogen	The antiserum was produced against synthesized peptide derived from human ZNF785. AA range:91-140
Specificity	ZNF785 Monoclonal Antibody detects endogenous levels of ZNF785 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	zinc finger protein 785;ZNF785
Observed Band	87kD
Calculated Molecular Weight	87kD
Cell Pathway	Nucleus .
Tissue Specificity	B-cell,Brain,Peripheral blood monocyte,Testis,
Function	function:May be involved in transcriptional regulation.,sequence caution:The predicted gene has been split into 2 genes: ZNF688 and ZNF785., similarity:Belongs to the krueppel C2H2-type zinc-finger protein family., similarity:Contains 1 KRAB domain.,similarity:Contains 2 C2H2-type zinc fingers.,similarity:Contains 7 C2H2-type zinc fingers.,
Background	function:May be involved in transcriptional regulation.,sequence caution:The predicted gene has been split into 2 genes: ZNF688 and ZNF785., similarity:Belongs to the krueppel C2H2-type zinc-finger protein family., similarity:Contains 1 KRAB domain.,similarity:Contains 2 C2H2-type zinc fingers.,similarity:Contains 7 C2H2-type zinc fingers.,


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using ZNF785 Monoclonal Antibody